

QUICK START GUIDE

PHOTON GENIE



Get Started

Rotate boxes and glass wands all over the body including major nerve and gland centers

(Energize the Body's Natural Healing Potential)

1. Plug in 2 transmitter boxes (*flat packs*) into main box
2. Plug in 2 glass wands into main box (one white/one red)
3. Plug main box into wall outlet
4. Turn main switch on (see green lights on small boxes)
5. Turn left switch on for glass wands (see pulsing lights)
6. Use intensity dial to set comfort level for wands
6. Begin with 20 minutes and work up to at least 1 hour
7. Continue increasing length of usage time as is comfortable
8. Sleep/rest with boxes as frequently as possible (more is better)
9. Hydrate well with water before and after

HEAD • NECK • CHEST • HANDS • FEET • STOMACH • UNDERARMS • GROIN •

PHOTON GENIE



2 GLASS WANDS:

(Use with or without the boxes)

- Place on bare skin on a single body
- This completes a Photon Genie *energy circuit*
- **Carpet shocks** may occur if the wands are not both touching bare skin on the same body
- Begin with the intensity dial at its lowest setting
- Adjust dial so it is comfortable for each body area (Different body parts have different sensitivity levels)
- Ozone is produced with the combination of skin contact and energy pulsing (Encouraging blood oxygenating)
- The intensity dial is only for the wands
- Turn wands off with the left ON/OFF switch when not in use and the intensity dial on lowest setting

2 TRANSMITTER BOXES:

(Use with or without the wands)

- Place smooth side without screws against the body
- Photon Genie energy flows through clothing, cloth, and bedding
- Does not have to touch the skin
- Green lights indicates power is on
- Green lights alternate between solid and flickering
(This is completely NORMAL)

BACK • GLANDS • LEGS • ARMS • SINUSES • EARS • TEMPLES • HIPS • KNEES

HEAD • NECK • CHEST • HANDS • FEET • STOMACH • UNDERARMS • GROIN • LEGS



SKILLING INSTITUTE

877.777.4788 | info@edskilling.com